


Version : 2021



画册的图片及参数仅供参考,如有更改,请按实际为准。

The picture and parameter are only for reference, will be based on the actual usage if any changes




Foshan Akos Medical Instrument Co., Ltd

佛山市耐可士医疗器械有限公司

3F, Unit A, No.4 Zone B, HAO SCIENCE PARK, Guicheng Street,  
Nanhai district, Foshan, Guangdong, China

佛山市南海区瀚天科技城B区4号楼 (A座3楼)

 (86+ / 0+ ) 757-86230689 / 757-86230699

 info@akosdent.com

 www.akosdent.com

**AKOS**



## ▶ Product contents

**K** Low speed handpiece

**A** High speed handpiece

**R** Handpiece spare part

**F** Niti File System

**B** Burs

## ABOUT US

### ● 关于我们

FOSHAN AKOS MEDICAL INSTRUMENT CO.,LTD is a professional dental handpiece and accessories manufacturer.

Our main products:High and low speed Handpiece,contra angles with most ratio..etc, we have an experienced R&D team, which can provide good OEM, ODM services. All of our products are CE & ISO certified.

Nowadays, we have established long term partnership with our clients from all over the world, and our products have been approved by many distributors and manufacturers.

AKOS DENT will never stop our innovative step and will continue to improve our products and quality in the future, in order to support our client with full effort, and provide the best after sales service to our clients.

Look forward to cooperation with you together !


佛山市耐可士医疗器械有限公司是一家专业生产牙科器械的企业。我们拥有一支专业、诚信、富有经验技术的团队。

我们专业生产：高速手机低速手机，各种比例弯机等，有专门经验丰富的研发团队，能提供定制服务和生产支持，所有产品均通过国内外重要认证。

目前，我们已经与世界各地的客户建立了长期的合作关系，我们的产品已经得到了许多经销商和制造商的认可。

我们将永远不会停止创新的步伐，并将继续精进我们的产品和质量，以全力支持我们的客户，并为我们的客户提供最好的售后服务。

我们期待与您的合作！





## Main Markets



Domestic Market

国内市场

Southeast Asia

东南亚

Europe

欧洲

Mid East

中东

South America

南美

- 2015** Established as a dental handpiece accessories manufacturer  
成立牙科配件生产线
- 2016** Established a handpiece technician develop team  
成立一支牙科手机技术研发队伍
- 2017** Starting to produce dental handpiece, and the sales is started growing  
开始生产和销售牙科手机, 并且销量逐渐增长
- 2018** Move to new factory, established a team for exportation, and a new brand AKOS  
搬迁至新厂房, 建立外贸团队, 国外品牌为AKOS
- 2019** Successfully got approval of the CFDA, starting to apply CE & ISO certifications  
成功取得国内医疗器械认证, 并开始申请CE和ISO等国际认证
- 2020** Successfully got approval of CE & ISO certifications, quality has been greatly improved  
成功取得CE和ISO等国际认证, 质量极大的提升, 完成业绩翻倍的目标

# LOW SPEED HANDPIECE



**1:5 / 1:4.2  
CONTRA ANGLE**



1:5 Contra Angle



German Bearing



Anti-Retraction System



Powerful & Quiet



Quatto Spray




# 1:4.2 CONTRA ANGLE

## DLC COATING & IMPORTED CERAMIC BALL BEARING



Gear and other spare parts are 100% compatible with the original Ti-Max Z45L, with the high quality DLC coating can ensure the long service life of our instruments, also without any uncomfortable noise to your patient. Imported ceramic bearing ensure the stability of our rotor, even under the extreme high speed and large load.



### Anti-retraction head, avoid cross-infection

- Special designed head, blood and dirt will face many obstacle to get inside the head, ensure the service life of the bearing and spindle



### Quattro spray, excellent cooling effect

- 4 points water spray, anti-retraction system, excellent cooling performance, high transparent fiber optic, ensure the light intensity



### DURA COAT

- Dura coat creates a smooth, scratch resistant surface which is comfortable to hold

## TECHNICAL PARAMETER 技术资料

Working dB 工作噪音	≤68dB
Working temperature 机身温度	≤20°C
Maximum speed 最高转速	Max. 168,000RPM
Head size 机头大小	w10mm/ h14.65mm




### High Accessibility

Unlike the common contra angle, this new design 45° surgical handpiece has a very high accessibility for the surgery, especially for the wisdom tooth.

## 1:4.2 CONTRA ANGLE

[1:4.2 增速弯机]



Reg. No. 注册号: KX95L (CFDA) / X95L (CE)

order-code 编号: **K25-45-CA-L** Fiber optic 光纤 P 按压式 **FG 高速 Ø 1.6**

## 1:4.2 EXTERNAL CONTRA ANGLE

[1:4.2 外水道增速弯机]



Reg. No. 注册号: KX95L (CFDA) / X95L (CE)


order-code 编号: **K25-45-CA-XL** Fiber optic 光纤 P 按压式 **FG 高速 Ø 1.6**

# 1:5 MINI HEAD CONTRA ANGLE

## DLC COATING & IMPORTED CERAMIC BALL BEARING



Gear and other spare parts are 100% compatible with the original Ti-Max X95, with the high quality DLC coating can ensure the long service life of our instruments, also without any uncomfortable noise to your patient. Imported ceramic bearing ensure the stability of our rotor, even under the extreme high speed and large load.



### Anti-retraction head, avoid cross-infection

- Special designed head, blood and dirt will face many obstacle to get inside the head, ensure the service life of the bearing and spindle



### Quattro spray, excellent cooling effect

- 4 points water spray, anti-retraction system, excellent cooling performance, high transparent fiber optic, ensure the light intensity



### DURA COAT

- Dura coat creates a smooth, scratch resistant surface which is comfortable to hold




### Special patent removable head design

- Special buckle design enable the head to be removable from its body without using any tool, easy for future maintenance

### Smaller and Slimmer

- This new model has a mini head and slim neck dimension compare with the old version and other competitors



### TECHNICIAN PARAMETER 技术资料

Working dB 工作噪音	≤68dB
Working temperature 机身温度	≤20°C
Maximum speed 最高转速	Max. 200,000RPM
Head size 机头大小	w8.8mm/ h12.1mm

## 1:5 MINI HEAD CONTRA ANGLE

[1:5 迷你增速弯机]



Reg. No. 注册号: KX95L (CFDA) / X95L (CE)

order-code 编号: **K25-CA-ML** Fiber optic 光纤 P 按压式 O FG 高速  $\phi$  1.6

# 1:5 INTERNAL / EXTERNAL CONTRA ANGLE

## TECHNICIAN PARAMETER 技术资料

Working dB 工作噪音	≤68dB
Working temperature 机身温度	≤20°C
Maximum speed 最高转速	Max. 200,000RPM
Head size 机头大小	w10mm/ h14.65mm

## 1:5 FULL S.S. EXTERNAL CONTRA ANGLE

[1:5 全不锈钢外水道增速弯机]



Reg. No. 注册号: KX95L (CFDA) / X95L (CE)

order-code 编号: **K25-CA-XL** **F** fiber optic 光纤 **P** 按压式 **O** FG 高速 **Ø 1.6**

## 1:5 FULL S.S. CONTRA ANGLE

[1:5 全不锈钢增速弯机]



Reg. No. 注册号: KX95L (CFDA) / X95L (CE)

order-code 编号: **K25-CA-L+** **F** fiber optic 光纤 **P** 按压式 **O** FG 高速 **Ø 1.6**

## 1:5 CONTRA ANGLE

[1:5 镀钛增速弯机]



Reg. No. 注册号: KX95L (CFDA) / X95L (CE), NX403 (CFDA) / X202 (CE)

order-code 编号: **K25-CA-L** **F** fiber optic 光纤 **P** 按压式 **O** FG 高速 **Ø 1.6**

order-code 编号: **K25-CA** **N** Non-fiber optic 非光纤 **P** 按压式 **O** FG 高速 **Ø 1.6**



# REDUCTION CONTRA ANGLE



20:1

20:1 Implant  
contra angle



German bearing



Powerful & Quiet



Quatto spray

75

Up to 75Ncm



# 20: 1 Implant CONTRA ANGLE

## TECHNICIAN PARAMETER 技术资料

Working dB 工作噪音	≤48dB
Working temperature 机身温度	≤20°C
Maximum speed 最高转速	2,000RPM
Head size 机头大小	w9.7 mm/ h13.78 mm

## 20:1 CONTRA ANGLE

[20:1 种植弯机]



Reg. No. 注册号: KX95L (CFDA) / X95L (CE), NX403 (CFDA) / X202 (CE)

order-code 编号: **K26-CA-L** **Fiber optic** 光纤 **P 按压式** **O CA减速 Ø 2.35**

order-code 编号: **K26-CA** **Non-fiber optic** 非光纤 **P 按压式** **O CA减速 Ø 2.35**

## > The handpiece can be used in any brand of implant system



### Large Torque

• Maximum torque required for operation is 65N, our designed maximum torque is 75N, the Minimum torque to break our handpiece is 125N



### Special designed irrigation

• The irrigation point is implanted inside the head, spray to the drill accurately, and not easy to be damaged while fall down to the floor





## > Developed for precision works



■ Durable gear box, can provide maximum 75Ncm for the surgery, it has a very strong transmission to provide maximum speed of 2000RPM, transmission using bearing instead of ball, performance is more smooth



■ The small head provide a better view for the doctor, making the surgery more comfortable and precise

## 20:1 LED GENERATOR CONTRA ANGLE

[20:1 风光种植弯机]



Reg. No. 注册号: KX95L (CFDA) / X95L (CE)

order-code 编号: **K26-CA-E** Generator 风光 P 按压式 O CA减速  $\varnothing$  2.35

## > Detachable body



■ Detachable bod, easy for clean after surgery and future maintainence

All the accessories are compatible with the original W&H contra angle WS-75LG

E-generator LED light, no need LED system in your implant motor

Powerful torque, can provide 75Ncm for surgery

## REDUCTION CONTRA ANGLE

### TECHNICIAN PARAMETER 技术资料

Working dB 工作噪音  $\leq 58$ dB

Working temperature 机身温度  $\leq 20^{\circ}\text{C}$

Gear Ratio 减速比 **4:1, 10:1, 16:1, 64:1**

## 4:1 CONTRA ANGLE

[4:1 弯机]



Reg. No. 注册号: KX95L (CFDA) / X95L (CE), NX403 (CFDA) / X202 (CE)

order-code 编号: **K27-CA-L** Fiber optic 光纤 P 按压式 O CA减速  $\varnothing$  2.35

order-code 编号: **K27-CA** Non-fiber optic 非光纤 P 按压式 O CA减速  $\varnothing$  2.35

## 16:1 CONTRA ANGLE

[16:1 弯机]



Reg. No. 注册号: KX95L (CFDA) / X95L (CE), NX403 (CFDA) / X202 (CE)

order-code 编号: **K30-CA-L** Fiber optic 光纤 P 按压式 O CA减速  $\varnothing$  2.35

order-code 编号: **K30-CA** Non-fiber optic 非光纤 P 按压式 O CA减速  $\varnothing$  2.35

## 4:1 CONTRA ANGLE

[4:1 弯机]



Reg. No. 注册号: KX95L (CFDA) / X95L (CE), NX403 (CFDA) / X202 (CE)

order-code 编号: **K27-CA-KL** **Fiber optic 光纤** **P 按压式** **CA 减速  $\varnothing$  2.35**

order-code 编号: **K27-CA-K** **Non-fiber optic 非光纤** **P 按压式** **CA 减速  $\varnothing$  2.35**

## 4:1 CONTRA ANGLE

[4:1 弯机]



Reg. No. 注册号: NK403 (CFDA) / N202 (CE)

order-code 编号: **K20-CA-4** **Non-fiber optic 非光纤** **P 按压式** **CA 减速  $\varnothing$  2.35**

## 4:1 POLISH CONTRA ANGLE

[4:1 抛光弯机]



Reg. No. 注册号: NK403 (CFDA) / N202 (CE)


order-code 编号: **K20-CA-4** **Non-fiber optic 非光纤** **P 按压式** **CA 减速  $\varnothing$  2.35**

## Endo Contra Angle

根管治疗用减速弯机

■ Fit to hand use files  
适用于手用根管锉

■ Reciprocate movement, imitate manual operation  
往复运动, 模仿手动操作



Hand use file  
Rotation speed : max 4,000rpm  
适用手用锉 / 转速最大达4,000rpm

## 10:1



Reg. No. 注册号: NK403 (CFDA) / N202 (CE)

order-code 编号: **K29-CA** **Non-fiber optic 非光纤** **P 按压式** **CA 减速  $\varnothing$  2.35**

## 64:1



Reg. No. 注册号: NK403 (CFDA) / N202 (CE)

order-code 编号: **K32-CA** **Non-fiber optic 非光纤** **P 按压式** **CA 减速  $\varnothing$  2.35**

## 16:1




Reg. No. 注册号: NK403 (CFDA) / N202 (CE)

order-code 编号: **K30-CA-H**

**Non-fiber optic 非光纤** **P 按压式** **CA 减速  $\varnothing$  2.35**

## 1:1



Reg. No. 注册号: NK403 (CFDA) / N202 (CE)

order-code 编号: **K20-CA-H**

**Non-fiber optic 非光纤** **P 按压式** **CA 减速  $\varnothing$  2.35**



13

135°C  
135°C high temperature  
& high pressure autoclavable

Can be washed with hot water

1 YEAR  
1 year warranty

135°C  
135°C high temperature  
& high pressure autoclavable

Can be washed with hot water

1 YEAR  
1 year warranty



14

# 4:1 Interproximal Enamel Kit

[4:1 切片弯机]



up and down movement  
上下运动

max10,000rpm  
Rotation speed 转速



- Pressure is more even, more accurate  
施压均匀, 更加精确
- Can be sterilized with high pressure  
可消毒 / 能承受高压

Reg. No. 注册号: NK403 (CFDA) / N202 (CE)

order-code 编号: **K27-CA-1** Non-fiber optic 非光纤 **CA减速  $\varnothing$  2.35**



Interproximal enamel kit  
[邻面去釉砂片套装]

Reg. No. 注册号: NK403 (CFDA) / N202 (CE)

order-code 编号: **K27-SET**



4:1 Reduction contra angle  
4:1 减速弯机



Gap measuring ruler  
测量尺  
0.10mm 0.15mm 0.20mm  
0.25mm 0.30mm 0.40mm  
0.50mm



Carborundum blade  
精细抛光刀片



Push button for remove  
卸针器



Manual Type Handle  
手动手柄



Milling strip  
砂条

Two side diamond (with hole) 有孔式(双侧切砂片)	15 $\mu$ m	15HD	0.17mm	2pcs	预抛光
	25 $\mu$ m	25HD	0.19mm	2pcs	精修
	40 $\mu$ m	40HD	0.23mm	2pcs	仿型
	60 $\mu$ m	60HD	0.32mm	2pcs	减速去釉
	90 $\mu$ m	90HD	0.40mm	2pcs	邻面开齿

1. Carborundum blade is stable and durable
2. It has a good compatibility with all kinds of handpiece
3. There are holes in the blade, which can effectively discharge the water and powder produced in the glaze.
1. 邻面片切砂片, 拥有金刚砂刀身, 持久耐用。
2. 采用了更加牢固的锁扣式设计, 保持完全固定和良好的兼容性。
3. 砂片上有孔, 可以有效排出去釉时产生的水和粉末。

# 1:1 HANDPIECE

CONTRA ANGLE /  
STRAIGHT HEAD /  
AIR MOTOR



# 1:1 Inner Channel CONTRA ANGLE

[1:1 内水道弯机]

## TECHNICIAN PARAMETER 技术资料

Working noise 工作噪音 ≤63dB Maximum speed 最高转速 200,000r/min

Air consumption 耗气量 56NL/min Air pressure 气压 2.8~3.0Kpa



Reg. No. 注册号: KX95L (CFDA) / X95L (CE)

order-code 编号: **K24-CA-G** Fiber optic 光纤 P 按压式 O FA 高速 Ø 1.6

order-code 编号: **K24-CA** Fiber optic 光纤 P 按压式 O CA 减速 Ø 2.35



Reg. No. 注册号: NX403 (CFDA) / X202 (CE)

order-code 编号: **K23-CA-G** Non-fiber optic 非光纤 P 按压式 O FA 高速 Ø 1.6

order-code 编号: **K23-CA** Non-fiber optic 非光纤 P 按压式 O CA 减速 Ø 2.35



Reg. No. 注册号: KX95L (CFDA) / X95L (CE), NX403 (CFDA) / X202 (CE)

order-code 编号: **K24-CA-K** Fiber optic 光纤 P 按压式 O CA 减速 Ø 2.35

order-code 编号: **K23-CA-K** Non-fiber optic 非光纤 P 按压式 O CA 减速 Ø 2.35

# 1:1 External CONTRA ANGLE

[1:1 外水道弯机]


- The accessories are compatible with the original NSK-Fx series contra angle FX25m.
- Suitable for the NSK-Fx Series air motor Fx205m (M4/B2)
- Durable stainless steel body

- 配件与原装NSK-Fx系列FX25m同型号弯机通配
- 适配NSK-Fx系列FX205m型号M4/B2马达
- 一体式不锈钢机头, 使用寿命长



Reg. No. 注册号: NK403 (CFDA) / N202 (CE)

order-code 编号: **K20-CA** Non-fiber optic 非光纤 P 按压式 O CA 减速 Ø 2.35



order-code 编号: > **K21-H**



compatible with/适配


**K21-CA**



Reg. No. 注册号: NK403 (CFDA) / N202 (CE)

order-code 编号: **K21-CA-P** Non-fiber optic 非光纤 P 按压式 O CA 减速 Ø 2.35

order-code 编号: **K21-CA** Non-fiber optic 非光纤 W 锁针式 O CA 减速 Ø 2.35



order-code 编号: > **K22-HP**



compatible with/适配


**K22-CA-P**



Reg. No. 注册号: NK403 (CFDA) / N202 (CE)

order-code 编号: **K22-CA-P** Non-fiber optic 非光纤 P 按压式 O CA 减速 Ø 2.35

order-code 编号: **K22-CA** Non-fiber optic 非光纤 W 锁针式 O CA 减速 Ø 2.35





## 1:1 Inner Channel/External STRAIGHT HEAD

[1:1 内/外水道直机]

### TECHNICIAN PARAMETER 技术资料

Working noise 工作噪音  $\leq 63\text{dB}$  Maximum speed 最高转速 200,000r/min

Air consumption 耗气量 56NL/min Air pressure 气压 2.8~3.0Kpa

### Inner Channel/内水道



Reg. No. 注册号: KX95L (CFDA) / K95L (CE), NK403 (CFDA) / N201 (CE)

order-code 编号: **K24-SH** **Fiber optic 光纤** **CA 减速  $\varnothing 2.35$**

order-code 编号: **K23-SH** **Non-Fiber optic 非光纤** **CA 减速  $\varnothing 2.35$**

### External/外水道



Reg. No. 注册号: KX95L (CFDA) / K95L (CE)

order-code 编号: **K24-SH-1** **Fiber optic 光纤** **CA 减速  $\varnothing 2.35$**



Reg. No. 注册号: NK403 (CFDA) / K201 (CE)

order-code 编号: **K22-SH** **Fiber optic 光纤** **CA 减速  $\varnothing 2.35$**

## 1:1 Inner/External STRAIGHT HEAD

[1:1 内/外水道直机]



Reg. No. 注册号: NK403 (CFDA) / K201 (CE)

order-code 编号: **K20-SH** **Non-Fiber optic 非光纤** **CA 减速  $\varnothing 2.35$**



Reg. No. 注册号: NK403 (CFDA) / K201 (CE)

order-code 编号: **K21-SH** **Non-Fiber optic 非光纤** **CA 减速  $\varnothing 2.35$**



Reg. No. 注册号: NK403 (CFDA) / K201 (CE)

order-code 编号: **K20-SH-1** **Non-Fiber optic 非光纤** **CA 减速  $\varnothing 2.35$**



Reg. No. 注册号: NK403 (CFDA) / K201 (CE)

order-code 编号: **K20-SH-2** **Non-Fiber optic 非光纤** **CA 减速  $\varnothing 2.35$**

## 1:1 Air Motor [1:1 气动马达]

### TECHNICIAN PARAMETER 技术资料

Working noise 工作噪音	≤63dB	Maximum speed 最高转速	20,000r/min
Air consumption 耗气量	56NL/min	Air pressure 气压	2.8~3.0 Kpa
Forward 正转	≥1.8000r/min	Reverse 反转	≥1.7000r/min



配置专利防堵过滤净水器,可调节水装置  
专利号:201620753014.5

### 气动马达特点: AIR MOTOR FEATURE:

- The whole body is made of stainless steel and imported EZO bearing. Self-lubrication eccentric stator, imported material molding blades
- Internal circulation cooling system, no vibration, no heat
- Patented water filter can effectively prevent the problems from the water spray and cooling system
- Compatible with all the handpiece that existing in the market

- 整体机身采用不锈钢外壳,进口EZO轴承自润滑偏心定子,进口材料成型叶片。
- 内循环冷却系统,不震动,不发热。搭配专利水过滤器结构,有效防止内水道水路堵塞等
- 冷却的不良问题。
- 通配市场所有低速内水道直机/弯机无缝对接。

## Inner Channel / 内水道



Reg. No. 注册号: KX95L (CFDA) / NX403 (CE)

order-code 编号: **K24-AM** **F**iber optic **光纤**



Reg. No. 注册号: NX203 (CFDA) / NX203 (CE), NX403 (CFDA) / NX403 (CE)

order-code 编号: **K23-AM** **N**on-fiber optic **非光纤**




## External/ 外水道




Reg. No. 注册号: NK203 (CFDA) / NK203 (CE), NK403 (CFDA) / NK403 (CE)

order-code 编号: **K22-AM** **N**on-fiber optic **非光纤**



Reg. No. 注册号: NK203 (CFDA) / NK203 (CE), NK403 (CFDA) / NK403 (CE)

order-code 编号: **K20-AM** **N**on-fiber optic **非光纤**



## 1:1 Inner/External LOW SPEED SET [1:1 内/外水道低速套装]

### Inner Channel/ 内水道

- The whole body is made of stainless steel, the head is manufactured by integral moulding, more durable.
- The body is according to human body mechanics, ultra silent, super torque.
- Patented water filter can effectively prevent the problems from the water spray and cooling system.

- 整体机身采用不锈钢材质,机头一体成型工艺,精工细琢、完美配备。
- 人体力学曲面外形,超静音超扭力。
- 内循环冷却,搭配专利水过滤器结构,有效防止内水道水路堵塞等冷却不良问题。



Reg. No. 注册号: KX95L (CFDA)

order-code 编号: **K24-SET** **F**iber optic **光纤** **P**按压式 **O**CA减速 **Ø 2.35**



# HIGH SPEED HANDPIECE



Reg. No. 注册号: NX203 (CFDA), NX403 (CFDA)

order-code 编号: **K23-SET** Non-fiber optic 非光纤 P 按压式 CA 减速  $\varnothing$  2.35



## External/ 外水道



Reg. No. 注册号: NK203 (CFDA), NK403 (CFDA)

order-code 编号: **K22-SET** Non-fiber optic 非光纤 W 锁针式 CA 减速  $\varnothing$  2.35



Reg. No. 注册号: NK203 (CFDA), NK403 (CFDA)

order-code 编号: **K21-SET** Non-fiber optic 非光纤 W 锁针式 CA 减速  $\varnothing$  2.35



Reg. No. 注册号: NK203 (CFDA), NK403 (CFDA)

order-code 编号: **K20-SET** Non-fiber optic 非光纤 P 按压式 CA 减速  $\varnothing$  2.35



# Fiber Optic Hanpiece

Compatible with NSK & KAVO



# Fiber Optic Hanpiece [高速光纤手机]

Compatible with NSK Machlite 通配NSK款光纤快接头

TECHNICIAN PARAMETER 技术资料	
Working dB 工作噪音	≤68dB
Water spraying mode 喷水方式	4 Points spray(四点喷水)
Maximum speed 最高转速	380,000~420,000rpm
Head size 机头大小	Standard (标准) w11.25/h13.3mm Torque (大头) w12.5/h13.3mm
	Mini (迷你) w10.08/h10.85mm Mini 45° (45°) w10.08/h10.85mm





Reg. No. 注册号: K9 (CFDA, CE)

order-code 编号: **Standard** **A06-SPQN**



Reg. No. 注册号: K9 (CFDA, CE)

order-code 编号: **Torque** **A06-TPQN**







Reg. No. 注册号: K9 (CFDA, CE)

order-code 编号: **Mini** **A06-MPQK**



Reg. No. 注册号: K9 (CFDA, CE)

order-code 编号: **Mini45°** **A06-45-MPQK**



### Excellent cooling performance

- 4 Point water & mist spray, anti-retraction, providing gentle mist during the surgery



### Imported ceramic ball bearing

- German steel plus japanese ceramic ball, combine a finest quality bearing for the turbine, ensure a stable performance and double the life cycle




### Compatibility

- Compatible with NSK type quick coupling, convenient for replacing the turbines




### Special patented filter system

- Equipped with a special filter design and one way valve, blocking most of the dirt from getting inside the water tube



### Powerful handpiece rotor

- High performance handpiece rotor, quiet and powerful with 18W, torque head with 25W



### High transparent fiber glass rod

- 40,000 LUX solid fiber optic cellular glass rod, ensure the sufficient illumination for the treatment


### NSK Type quick coupling



# Fiber Optic Hanpiece [高速光纤手机]


Compatible with KAVO Multiflex 通配KAVO款光纤快接头

TECHNICIAN PARAMETER 技术资料	
Working dB 工作噪音	≤68dB
Water spraying mode 喷水方式	4 Points spray(四点喷水)
Maximum speed 最高转速	380,000~420,000rpm
Head size 机头大小	Standard (标准) w11.25/h13.3mm Torque (大头) w12.5/h13.3mm
	Mini (迷你) w10.08/h10.85mm Mini 45° (45°) w10.08/h10.85mm



Reg. No. 注册号: K9 (CFDA, CE)

order-code 编号: **Standard** **A06-SPQK**



Reg. No. 注册号: K9 (CFDA, CE)

order-code 编号: **Torque** **A06-TPQK**



Reg. No. 注册号: K9 (CFDA, CE)

order-code 编号: **Mini** **A06-MPQK**



Reg. No. 注册号: K9 (CFDA, CE)

order-code 编号: **Mini 45°** **A06-45-MPQK**

## Excellent cooling performance

- 4 Point water & mist spray, anti-retraction, providing gentle mist during the surgery

## Imported ceramic ball bearing

- German steel plus japanese ceramic ball, combine a finest quality bearing for the turbine, ensure a stable performance and double the life cycle



# LED Handpiece with Generator



## Compatibility

- Compatible with KaVo Multiflex quick coupling, convenient for replacing the turbines




## Special patented filter system

- Equipped with a special filter design and one way valve, blocking most of the dirt from getting inside the water tube



## Powerful handpiece rotor

- High performance handpiece rotor, quiet and powerful with 18W, torque head with 25W



## High transparent fiber glass rod

- 40 · 000 LUX solid fiber optic cellular glass rod, ensure the sufficient illumination for the treatment

## KAVO Type quick coupling



order-code 编号:



**Q-K2-M6**



6 hole 六孔



order-code 编号:



**Q-K3-M4**



4 hole 自发光 / 四孔 **Generator**



order-code 编号:



**Q-K3-B2**



2 hole 自发光 / 二孔 **Generator**







order-code 编号:  
**R05-MP**  
compatible with/适配  
**A05-MPQ**



**Couplings**

Reg. No. 注册号: X7 (CFDA) / K9 (CE)

**2 LED Paediatric 双灯**

order-code 编号: **Mini** | **A05-MPQ** | 1 way spray 4 air | 单点喷水



order-code 编号:  
**R05-45SP**  
compatible with/适配  
**A05-45-MP**



**Couplings**

Reg. No. 注册号: X7 (CFDA) / K9 (CE)

**Single LED 单灯**

order-code 编号: **Mini45°** | **A05-45-MP** | 1 way spray | 单点喷水



order-code 编号:  
**R05-SP**  
compatible with/适配  
**A05-SP**



**Couplings**

Reg. No. 注册号: X7 (CFDA) / K9 (CE)

**Single LED 单灯**

order-code 编号: **Standard** | **A05-SP** | 1 way spray 1 air | 单点喷水



order-code 编号:  
**R05-TP**  
compatible with/适配  
**A05-TP**



**Couplings**

Reg. No. 注册号: X7 (CFDA) / K9 (CE)

**Single LED 单灯**

order-code 编号: **Torque** | **A05-TP** | 1 way spray | 单点喷水



### Powerful Torque Cartridge

- With imported bearing, rotor work more stable, silence, and life cycle is much longer



### Finest craft spindle

- We control all concentricity less than 15µm, some can be less than 10µm, we don't have problem of losing bur



### Powerful generator

- Generator is the heart of the LED, our generator can provide a high illumination, and at the same time more durable against the high temperature sterilization



### Compatible with the quick coupling


- Our generator handpiece can work with the international standard coupling, more convenient during the surgery

# Standard Handpiece



## 4 Points spray handpiece [4点喷水手机]

order-code 编号:  
**R04-SP**  
compatible with/适配  
**A04-SP**



OO or OO  
2 hole 4 hole

Reg. No. 注册号: K1 (CFDA, CE), K6 (CFDA, CE)

order-code 编号: **Standard** **A04-SP**



order-code 编号:  
**R04-TP**  
compatible with/适配  
**A04-TP**




OO or OO  
2 hole 4 hole

Reg. No. 注册号: K1 (CFDA, CE), K6 (CFDA, CE)

order-code 编号: **Torque** **A04-TP**





order-code 编号:  
**R04-MP**  
compatible with/适配  
**A04-MP**



OO or OO  
2 hole 4 hole

Reg. No. 注册号: K1 (CFDA, CE), K6 (CFDA, CE)

order-code 编号: **Mini** **A04-MP**



order-code 编号:  
**R04-45-SP**  
compatible with/适配  
**A04-45-SP**



OO or OO  
2 hole 4 hole

Reg. No. 注册号: N5 (CFDA, CE)

order-code 编号: **Mini45°** **A04-45-SP**



order-code 编号:  
**R04-SP**  
compatible with/适配  
**A04-SPQ**



OO or OO  
2 hole 4 hole

Reg. No. 注册号: N8 (CFDA, CE)

order-code 编号: **Standard** **A04-SPQ**



**Couplings**



order-code 编号:  
**R04-TP**  
compatible with/适配  
**A04-TPQ**



OO or OO  
2 hole 4 hole

Reg. No. 注册号: N8 (CFDA, CE)

order-code 编号: **Torque** **A04-TPQ**



**Couplings**



## Anti-retraction handpiece 防回吸卫生机头手机

### Cartridge compatible with NSK

- Special design anti-retraction head, prevent the contamination from getting in the head
- Avoiding cross-infection
- Compatible with NSK® Pana Max
- 特殊卫生机头设计, 防止杂质进入机头
- 避免交叉感染
- 机芯适配于 NSK Pana Max



### TECHNICIAN PARAMETER 技术资料

Working dB 工作噪音	≤68dB
Power 功率	8.5N
Maximum speed 最高转速	350,000~380,000rpm
Head size 机头大小	Standard (标准) w11.25/h13.3mm Torque (大头) w12.5/h13.3mm

order-code 编号:  
**R02-SPM1**  
compatible with/适配  
**A02-SPM1**



OO or OO  
2 hole 4 hole

Reg. No. 注册号: K1 (CFDA, CE), K6 (CFDA, CE)

order-code 编号: **Standard** **A02-SPM1** | 1 way spray 3 air | 单点喷水3点喷气



Reg. No. 注册号: K1 (CFDA, CE), K6 (CFDA, CE)



order-code 编号: **Torque** **A02-TPM1** **1 way spray 3 air** | 单点喷水3点喷气



Reg. No. 注册号: K1 (CFDA, CE), K6 (CFDA, CE)



order-code 编号: **Standard** **A02-SPMQ1** **1 way spray 3 air** | 单点喷水3点喷气



Reg. No. 注册号: K1 (CFDA, CE), K6 (CFDA, CE)



order-code 编号: **Torque** **A02-TPMQ1** **1 way spray 3 air** | 单点喷水3点喷气

## Handpiece Kit

高低速手机组合套装/学生套装




### CUSTOMIZE YOUR OWN KIT

#### > AK31-KIT-X

##### 1 or 2 High 1 Low complete set

- High speed handpiece x 1 or 2
- Contra Angle x 1
- Straight head x 1
- Air motor x 1
- Luxury Aluminum or Normal Plastic

**Free** Cartridge for high speed x 1



#### > AK31-KIT-1

##### 1 High 1 Low complete set

- High speed handpiece (Push button) x 1
- Contra Angle (Latch type) x 1
- Straight head x 1
- Air motor (External water) x 1

**Free** Cartridge for high speed x 1



# HANDPIECE ACCESSORIES



## > AK31-KIT-2

### 2 High 1 Low complete set

- High speed handpiece (Wrench type) x 2
- Contra Angle (Latch type) x 1
- Straight head x 1
- Air motor (External water) x 1

**Free** Cartridge for high speed x 1



## > AK31-KIT-3

### 2 High 1 Low complete set

- High speed handpiece (Push button) x 2
- Contra Angle (Latch type) x 1
- Straight head x 1
- Air motor (External water) x 1

**Free** Cartridge for high speed x 1



## > AK31-KIT-4

### 1 High 1 Low complete set

- Anti-Retraction handpiece (Push button) x 1
- Contra Angle (Latch type) x 1
- Straight head x 1
- Air motor (External water) x 1

**Free** Cartridge for high speed x 1



## > AK31-KIT-5

### 1 High 1 Low complete set

- Quattro spray handpiece (Push button) x 1
- Contra Angle (Push button) x 1
- Straight head x 1
- Air motor (Inner water) x 1

**Free** Cartridge for high speed x 1



# Handpiece spare part

[牙科手机配件]

R

## Low speed part

### Shaft 传动轴

1:1

> Shaft 1-1

16:1

> Shaft 16-1

20:1

> Shaft 20-1

1:5

> Shaft 1-5

1:4.2

> Shaft 1-5

### Cartridge 机芯



> R 1-1

1:1



> R 16-1

16:1



> R 20-1

20:1



> R 1-5

1:5



> R 1-4.2

1:4.2

### Head cap of contra angle 后盖



> Cap 1-1

1:1



> Cap 16-1

16:1



> Cap 20-1

20:1



> Cap 1-5

1:5



> Cap 1-4.2

1:4.2

## High speed part

### Cartridge 机芯

FOR NSK®



NPA-M03

> R01-M



NPA-MU03

> R01-MP



NPA-S03

> R01-S



NPA-SU03

> R01-SP



NPA-T03

> R01-T



NPA-TU03

> R01-TP



PAX-SU03

> R02-SPM1



PAX-TU03

> R02-TPM1



PAX2-SU03

> R01-SPM2



SX-MU03

> R01-SPM3



SX-SU03

> R02-TPM3



NSK-X500L

> R-X500L



NSK-X600L

> R-X600L

FOR KaVo®



KAVO-8000

> R-8000



KAVO-105/660

> R-105



KAVO-637

> R-637



原版

KAVO-645/642 Original

> R-645 Original



巴西版

KAVO-636/646 Brasil

> R-636 Brasil



德国版

KAVO-636/640/625/655 Germany

> R-636 Germany

FOR Sirona®



标准版

SIRONA T3 standard

> R-T3standard



迷你版

SIRONA T3 mini

> R-T3mini

FOR W&H®



W&H TE-96

> R-TE96



W&H TE-98






> R-TE98






## Shaft of cartridge high speed handpiece 高速手机机芯按压轴

Push button spindle 按压款

Wrench button spindle 锁针款

-  10.6 mm [PS-10.6](#)
-  11.8 mm [PS-11.8](#)
-  12.6 mm [PS-12.6](#)
-  12.75 mm [PS-12.75](#)
-  13.75 mm [PS-13.75](#)

-  10.85 mm [WS-10.85](#)
-  12.55 mm [WS-12.55](#)
-  13.25 mm [WS-13.25](#)



 [FIT IN PAX-SU03](#)

order-code  
**Group-PAX-SU03**

 [FIT IN PAX-S03](#)

order-code  
**Group-PAX-S03**



 [FIT IN PAX-TU03](#)

order-code  
**Group-PAX-TU03**

 [FIT IN PAX-T03](#)

order-code  
**Group-PAX-T03**




 [FIT IN PAX2-SU03](#)

order-code  
**Group-PAX2-SU03**

 [FIT IN PAX2-S03](#)

order-code  
**Group-PAX2-S03**




 [FIT IN SX-MU03](#)

order-code  
**Group-SX-MU03**



 [FIT IN SX-SU03](#)

order-code  
**Group-SX-SU03**



 [FIT IN NPA-MU03](#)

order-code  
**Group-NPA-MU03**

 [FIT IN NPA-M03](#)

order-code  
**Group-NPA-M03**



 [FIT IN NPA-SU03](#)

order-code  
**Group-NPA-SU03**

 [FIT IN NPA-S03](#)

order-code  
**Group-NPA-S03**



 [FIT IN NPA-TU03](#)

order-code  
**Group-NPA-TU03**



 [FIT IN NPA-T03](#)

order-code  
**Group-NPA-T03**

## Head cap of high speed handpiece 高速手机后盖

Fit into NSK®



Fit in Sirona®



Fit into KaVo®

## Quick coupling of high speed handpiece 高速手机快接头

Fit into KaVo®




Fit into NSK®









## Impeller specification 风轮

Fit into NSK®




## Fit into W&H®




## Bearing of high speed handpiece 轴承



## Head spare part of speed handpiece 低速手机机头组件



## Cartridge repair kit 手机维修工具



# Niti File System

[镍钛锉系统]

F

## Niti Protaper (Blue)

镍钛大锥度蓝色锉

### > F-001

- CM-Niti material with heat treatment technology, more flexible and durable
- Security tip, reduces the possibility of root canal deviation
- Smart varied taper design, high efficiency
- Made of high quality controlled memory wire, after heat activation treatment blue color is derived
- Blue file is flexible, tough, and pre-curved suitable for curved canals
- More sharper and performs better for canal preparation
- 6pcs per pack, available for single or assorted size



Controlled Memory  
Thermal Activation  
Triangle Cross Section  
Anti-Broken

Product Name 产品名称	Length 长度(mm)	Size 型号	Cross-Sectional 横截面	Taper 锥度 (%)	Tip 针尖尺寸(mm)	Recommended rotary speed 推荐转速	Torque 扭力	Packing 包装规格				
Niti Protaper (Blue) 镍钛大锥度蓝色锉	19	○ SX		3.5-9	0.19	150-300 RPM	2.0-2.5 N.cm	6pcs/box or blister Single or Assorted are available 可选单/混装				
		● S1		2-11	0.17							
	● S2	4-11.5		0.20								
	● F1	7		0.20								
	● F2	8		0.25								
	● F3	9		0.30								
	● F4	6		0.40								
	● F5	5		0.50								
	21/25/31											

## Preshaper

D3 锉

### > F-002

- High quality of Niti Alloy, after special heat treatment, the file is more flexible, almost impossible to break
- Special design has many advantages for cutting, and easy to remove the debris
- Safe round shape of the top, avoid creating step or side perforation
- Rotating smooth and responsively in the canal, more efficient and strong
- 6pcs per pack, available for single or assorted sizes  
(2pcs hand K files, 4pcs engine files)



Product Name 产品名称	Length 长度(mm)	Size 型号	Recommended rotary speed 推荐转速	Torque 扭力	Packing 包装规格
Preshaper D3锉	21/25/31	● S1	250-300 RPM	2.5N.cm	6pcs/box Single or Assorted are available 可选单/混装
		● S2			
		● S3			
		R	350 RPM	2.5N.cm	

## Niti Protaper (Gold, Thermal Activation)

镍钛机用锉 (金色热激活)

### > F-003

- The file is multi-tapered, more flexible and has a better performance while cutting
- High quality Niti Alloy, will become gold color after heat treatment
- 6pcs per pack, available for single or assorted sizes.



Advanced heat treatment  
Special processing technique  
More Sharp Sharper  
Strong & Flexible

Product Name 产品名称	Length 长度(mm)	Size 型号	Cross-Sectional 横截面	Taper 锥度 (%)	Tip 针尖尺寸(mm)	Recommended rotary speed 推荐转速	Torque 扭力	Packing 包装规格				
Niti Protaper (Gold, Thermal Activation) 镍钛机用锉 (金色热激活)	19	○ SX		3.5-9	0.19	150-300 RPM	2.0-2.5 N.cm	6pcs/box or blister Single or Assorted are available 可选单/混装				
		● S1		2-11	0.17							
	● S2	4-11.5		0.20								
	● F1	7		0.20								
	● F2	8		0.25								
	● F3	9		0.30								
	● F4	6		0.40								
	● F5	5		0.50								
	21/25/31											

## Kids File (for children)

乳牙锉

### > F-004

- Short thread design, matching the short root canal of baby's teeth;
- The small files protects the fragile root canal wall of baby's teeth effectively;
- The thread is changed from the tip of the files to the back to facilitate chip evacuation




Product Name 产品名称	Length 长度(mm)	Size 型号	Cross-Sectional 横截面	Taper 锥度 (%)	Tip 针尖尺寸(mm)	Recommended rotary speed 推荐转速	Torque 扭力	Packing 包装规格
Kids Files 乳牙锉	12mm	25#04		4	0.25	150-300 RPM	2.0-2.5 N.cm	6pcs/box Single or Assorted are available 可选单/混装
	12mm	25#06		6	0.25			
	12mm	30#04		4	0.25			
	14mm	40#04		4	0.25			

### SUPER FLEXI FILES (Thermal Activation)

超韧锉 (热激活功能)

> **F-005**

- PRO-FLEXI FILES is a simple but effective rotary system.
- Special heat treatment result in blue color
- One sequence for all root canals




Product Name 产品名称	Length 长度(mm)	Size 型号	Cross-Sectional 横截面	Taper 锥度 (%)	Tip 针尖尺寸(mm)	Packing 包装规格
<b>SUPER FLEXI FILES</b> (Thermal Activation) 超韧锉 热激活功能	19	01/09		9	0.25	6pcs/box or blister Single or Assorted are available 可选单/混装
		20/04		4	0.20	
		25/04		4	0.25	
		30/04		4	0.30	
		35/04		4	0.35	
		25/06		6	0.25	

### One File (compatible with Reciproc VDW)

单支锉

> **F-006**

- Only one file can finish preparing a root canal
- With less procedure, effective and time saving
- Special thread design, reciprocating movement
- 3 different sizes suitable for all kinds of root canals
- 6pcs per pack, available for single or assorted sizes.



Product Name 产品名称	Length 长度(mm)	Size 型号	Cross-Sectional 横截面	Taper 锥度 (%)	Tip 针尖尺寸(mm)	Packing 包装规格
<b>One File</b> (compatible with Reciproc VDW) 单支锉	21/25/31	R25		8	0.25	6pcs/box or blister Single or Assorted are available 可选单/混装
		R40		6	0.40	
		R50		5	0.50	

## FG diamond burs [金刚砂车针]

B



### Taper Fissure

Crown and bridge preparation



	3.5 19	3.5 19	4.0 19	4.0 19	7.0 21	10 24	4.0 19	7.0 21
▲	TF-009	TF-011	TF-012	TF-014S	TF-014	TF-014EL	TF-016S	TF-016
JP Type	-	TF-41	TF-42	TF-43	TF-20	TF-11	TF-31	TF-21
EU Type	845009	845011	845012	845014	847014	848014	845016	847016

	10 22	5.0 20	10 22	5.5 19	10 22	4.0 19	4.0 19	4.0 19
▲	TF-016L	TF-018	TF-018L	TF-021	TF-023L	TF-016R	TF-018R	TF-016
JP Type	TF-12	TF-23	TF-13	TF-22	TF-14	TF-11	TF-31	TF-21
EU Type	848016	845018	848018	845021	848023	845KR016	845KR018	847016



### Wheel

Crown and bridge preparation  
Occlusal reduction



	1.8 19
▲	WR-042
JP Type	WR-13
EU Type	909042




### Conical Ogival

Crown and bridge preparation  
Cavity Preparation



	6 19	10 24	7.0 21	7.0 21
▲	CO-012	CO-012EL	CO-014	CO-016
JP Type	F0-41	F0-11	F0-21	F0-22
EU Type	861012	863012	862014	862016




### Torpedo

Cavity Preparation




	6.0 19	6.0 21	8.0 22	8.0 20	8.0 22	8.0 22
▲	TB-012S	TB-012	TB-014L	TB-016	TB-016L	TB-018L
JP Type	S0-S20	S0-20	S0-21	-	-	-
EU Type	-	845011	845012	-	847014	848014





### Round End taper fissure

Crown and bridge preparation



▲ TR-012	TR-014	TR-016S	TR-016	TR-016L	TR-A016S	TR-016
JPType	TR-20	TR-S21	TR-21	TR-11	TR-12	TR-13
EUType	856012	856014	856016	856016	850016	856018

▲ TF-021L	TF-021EL	TF-023EL	
JPType	TR-14	TR-15	TR-19
EUType	856021	850021	850023

### Taper Conical


Crown and bridge preparation  
Occlusal reduction / Cavity Preparation



▲ TR-012	TR-014	TR-016S	TR-016
JPType	TR-20	TR-S21	TR-21
EUType	856012	856014	856016

### Flat Fissure

Crown and bridge preparation



▲ FF-010	FF-012L	FF-013	FF-014	FF-014L	FF-016L
JPType	SF-11	SF-31	SF-21	SF-12	SF-13
EUType	835010	837012	835013	836014	837014

### Round

Cavity Preparation



▲ RB-006	RB-008	RB-010	RB-012	RB-014	RB-016
JPType	BR-49	BR-45	BR-46	BR-41	BR-40
EUType	801006	801008	801010	801012	801014

▲ RB-018	RB-023	
JPType	BR-31	-
EUType	801018	801023

### Pear


Cavity Preparation



▲ PR-008	PR-012	PR-018	PR-021	PR-032
JPType	-	EX-20	EX-21	EX-26
EUType	830008	830012	830018	830021

### BUD

Cavity Preparation



▲ BB-016	BB-018	BB-028	BB-032	EB-018	EB-023
JPType	-	F0-25	F0-27	EX-26	EX-40
EUType	390016	368018	368028	368032	379023



## FG Carbide burs

[FG钨钢车针]

B

**Round**

US No.	1/4	1/2	1	2	3	4
FG 500 314 001 006...	005	006	008	010	012	014
Order No.	RB-005	RB-006	RB-008	RB-010	RB-012	RB-014

**Round End Fissure**

US No.	1555	1556	1557	1558	1559
FG 500 314 137 008...	006	008	009	010	012
Order No.	REF-006	REF-008	REF-009	REF-010	REF-012

**Pear**

US No.	329	330	331	332	333
FG 500 314 237 006...	006	008	009	010	012
Order No.	PR-006	PR-008	PR-009	PR-010	PR-012

**Taper Fissure**

US No.	698	699	700	701	702
FG 500 314 168 008...	008	009	010	012	014
Order No.	TF-SB008	TF-SB009	TF-SB010	TF-SB012	TF-SB014

**Inverted Cone**

US No.	33 1/2	34	35	36	37
FG 500 314 010 006...	006	008	010	012	014
Order No.	INC-006	INC-008	INC-010	INC-012	INC-014

**Flat End Fissure**

US No.	557	558	559
FG 500 314 107 008...	010	012	014
Order No.	FF-B010	FF-B012	FF-B014

**Flat End Fissure**

US No.	55	57	58
FG 500 314 107 006...	008	010	012
Order No.	FF-S008	FF-S010	FF-S012

## HP Carbide burs

[HP低速直柄钨钢车针]

B

**Round**

US No.	1/4	1/2	1	2	3	4	5	6	7	8
HP 500 104 001 006...	005	006	008	010	012	014	016	018	021	023
Order No.	RB-005	RB-006	RB-008	RB-010	RB-012	RB-014	RB-016	RB-018	RB-021	RB-023

**Round End Fissure**

US No.	1555	1556	1557	1558	1559	1560	1561	1563
HP 500 104 137 008...	006	008	009	010	012	014	016	018
Order No.	REF-006	REF-008	REF-009	REF-010	REF-012	REF-014	REF-016	REF-018

**Pear**

US No.	329	330	331	332	333
HP 500 104 237 006...	006	008	009	010	012
Order No.	PR-006	PR-008	PR-009	PR-010	PR-012

**Taper Fissure**

US No.	698	699	700	701	702	703
HP 500 104 168 008...	008	009	010	012	014	016
Order No.	TF-SB008	TF-SB009	TF-SB010	TF-SB012	TF-SB014	TF-SB016

**Inverted Cone**

US No.	33 1/2	34	35	36	37	38	39	40	41
HP 500 104 010 006...	006	008	010	012	014	016	018	021	023
Order No.	INC-006	INC-008	INC-010	INC-012	INC-014	INC-016	INC-018	INC-021	INC-023

**Flat End Fissure**

US No.	557	558	559	560	561	562	563
HP 500 104 107 008...	010	012	014	016	018	021	023
Order No.	FF-B010	FF-B012	FF-B014	FF-B016	FF-B018	FF-B021	FF-B023

**Flat End Fissure**

US No.	55	57	58	60
HP 500 104 107 006...	008	010	012	016
Order No.	FF-S008	FF-S010	FF-S012	FF-S016



## OEM/ODM service ✓ 订制服务

### 1 ODM SERVICE (MOQ 300~500PCS)

1 Completely your own or compatible with any preferable brand

NSK® KaVo® W&H

2 Spray method  
\_\_\_\_\_ way spray



3 Pattern of the handpiece



Others

4 Material of the handpiece



Copper



Stainless Steel




Titanium

5 Private logo



### 2 OEM SERVICE (NO MOQ REQUIRED)

Your private Logo will be laser printed on our product



We can also supply the above products for you, if there's further you need, please feel free to E-mail us, we will provide all details including catalogue and price list to you!



# WARRANTY TERM

## 售后保障

We will provide warranty and after sales service for all our products, for those A+ quality product we will provide at least 1 Year Warranty.

我们会对我们的产品提供售后服务, A+质量的产品会提供至少1年的保修期。

We have after sales service center in Europa and China, if there' s any quality issue while using our product, as long as the products are still within the warranty period, and we will provide repair and maintenance free of charge.

However, any inappropriate operation that cause damage to the instrument, will not covered by the warranty.

我们在欧洲和中国都有售后维修中心, 如果使用过程中遇到任何质量问题, 只要产品仍然在保修期内, 我们都可以提供免费的保修服务。如果产品是由于不当的使用而造成损坏, 则不属于保修范围内。

To ensure the long service life of the instrument, please maintain the instrument in time, most of the defect problems are mainly due to the inappropriate maintenance, here are some common cases that result in failure of the instrument.

确保产品有更长的使用周期, 请及时对产品进行维护, 很大一部分产品损坏的情况都是由于不当的维护或者没有对产品进行维护而造成的, 下面是一些常见的容易导致产品损坏的实例

Instrument without maintenance



(没有保养的设备)

Incorrect or cheap lubrication oil (Dark and dirty color)



(不正确的或使用劣质的手机油)



# AKOS